

Manorapid® Synergy



Alcoholic solution for the hygienic and surgical hand disinfection for areas with increased risk of contamination with viruses.

Range of application:

There are viruses that can be inactivated comparably simple and there are viruses (i.e. the polio virus) whose destruction has been posing an almost insuperable conflict to microbiologists so far. Although high percentage alcohol (i.e. Ethanol 90–95%) offers the desired effect, it stresses the skin and can only be stored in small quantities due to its low flashpoint. Manorapid® Synergy is a completely new formulation, a new chemical efficacy mechanism, made possible by discovery of the effect that alcohols in low concentration can be significantly increased through additives, so called "synergists". Thus, a significantly rapid efficacy is possible against coated and uncoated viruses in max. 1 min. The formulation is already protected by patent in several countries. Manorapid® Synergy is a ready for use alcoholic solution for the hygienic and surgical hand disinfection.

Special properties:

- broad spectrum of efficacy
- strong caring characteristic
- pleasant scent

Use:

● Hygienic hand disinfection:

Manorapid® Synergy is already effective after 15 sec. Upon current recommendations we suggest a contact time of 30 sec.

● Virus inactivation:

Rub into hands undiluted and keep them moist for 1 min.

● Surgical hand disinfection:

Rub into hands undiluted and keep them moist for min. 1.5 min.

Spectrum of efficacy:

- Bactericidal (incl. M. terrae)
- Levurocidal
- Virucidal*/**

Active ingredients:

100 g contain:
10 g 1-Propanol, 57.6 g Ethanol (96%)

Listing:

VAH/ÖGHMP/RKI

Physical and chemical characteristics:

Appearance: clear, colorless liquid
Odour: alcoholic
pH Value: ca. 3 at 20°C
Density: 0.88 g/cm³ at 20°C
Flashpoint: 22°C DIN 51755

Presentation:

- Cartons with 30 x 150 ml, 20 x 500 ml and 10 x 1 L plastic bottles
- 5 L plastic containers
- Wall dispensers deliverable for 500 ml and 1 L bottles

* acc. to RKI-recommendation Bundesgesundheitsblatt 01/2004 (Polio, Adeno, SV40 and Vaccinia)
** EN 14476 (2005) (Polio and Adeno)

Hygienic hand disinfection (incl. M. terrae)	3 ml - 30 sec.
--	----------------

Surgical hand disinfection (incl. Mycobacteria)	2x 3ml - min. 1.5 min.
---	------------------------

Limited virucidal (incl. HIV, HBV, HCV)*	30 sec.
--	---------

Virucidal* (accord. to RKI)	1 min.
-----------------------------	--------

Virucidal** (accord. to EN 14476)	30 sec.
-----------------------------------	---------

SV40	1 min.
------	--------

Vaccinia	30 sec.
----------	---------

BVDV	30 sec.
------	---------

Rotavirus	15 sec.
-----------	---------

Influenza A and B	30 sec.
-------------------	---------

HSV1 and HSV2	30 sec.
---------------	---------

SARS	30 sec.
------	---------

HAV	30 sec.
-----	---------

FCV (surrogate virus for Noro virus)	15 sec.
--------------------------------------	---------

This product is manufactured according to the European Quality Standard ISO 9001. Further information upon request by your expert adviser (expert opinions) or on our home page.

Use biocides safely. Always read labelling and product information before use.

Made in Germany

ANTISEPTICA



ANTISEPTICA chem. pharm. Produkte GmbH
D-50259 Pulheim, Carl-Friedrich-Gauß-Straße 7
Tel.: +49-2234-984 66-0, Fax: +49-2234-984 66-11
E-Mail: info@antiseptica.com, www.antiseptica.com

QUALITÄTS-MANAGEMENT
Wir sind zertifiziert
Regelmäßige freiwillige
Überwachung nach EN ISO 13485:2003



QUALITÄTS-MANAGEMENT
Zertifiziert nach
Medizinprodukte-Richtlinie
93/42/EWG



QUALITÄTS-MANAGEMENT
Wir sind zertifiziert
Regelmäßige freiwillige
Überwachung nach ISO 9001:2000



ANTISEPTICA Austria GmbH
A-1020 Wien, Wehlstraße 150/Haussteinstraße, Top 61
Tel.: +43-1-374 66 00, Fax: +43-1-374 66 00-66
E-Mail: office@antiseptica.at, www.antiseptica.at

Antisepsis by Antiseptica – Antisepsis mit Vernunft